

For claw treatment:

Demotec FuturaPad

The revolution in claw treatment



Demotec FuturaPad

The revolution in claw treatment

Demotec has developed a unique system for claw treatment for which a patent application has been made – covering most major countries of the world. FuturaPad makes claw treatment more effective; it is simple fast and easy to use.

FuturaPad is unique; one surface of the wooden block employed is covered with a flexible reactive Pad, which overcomes any unevenness of the sole and maximises adhesion over the full surface area. Bearing in mind the anatomy and the function of the claw, the reactive Pad is shaped correctly and is affixed to the wooden block in the optimal position. FuturaPad's reactive Pad is activated by a special glue presented in a unit application pack.

The processing is done without any mixing; all components needed are included in the unit application pack and the rate of polymerisation is not dependent on the temperature; which is of considerable benefit, particularly in the winter.



- **Activate the Pad, press the block onto the claw - done!**
- **10-15 seconds later the block is fixed firmly onto the claw**
- **Polymerisation within seconds – even during the winter**



- **Presented in a unit application pack**
- **No extra components needed**
- **Available in two sizes (normal or XL)**



- **FuturaPad is shaped to conform to the anatomy of the claw**
- **The flexible Pad overcomes any unevenness of the sole**
- **No health transport or storage problems**

Demotec – The competent experts in claw treatments

Demotec Siegfried Demel

Brentanostr. 21

D – 61130 Nidderau

Germany

Telefon +49-6187-21200, Telefax +49-6187-21208

Internet www.demotec.com, Email demotec@demotec.com

Demotec[®]

Intelligent Problem Solving